

SDO[®] Lite

For other languages, please download from: www.jobskin.co.uk/wash

How to care for your new ready-to-wear garment

Before washing your garment, please follow these simple rules to help protect and prolong the life of your garment:

REMEMBER

- Wash your garment 2-3 times a week to rejuvenate the fabric
- Turn inside out and fasten zips and hook and loop tabs
- Wash dark colours separately
- Do not dry in direct sunlight or on a radiator
- **DO NOT USE A FABRIC CONDITIONER OR A DETERGENT THAT CONTAINS AN ADDED FABRIC CONDITIONER**

HAND WASH

- If hand washing, roll your garment inside a thick towel to remove any excess moisture
- Rinse well and dry flat on a clean, dry towel. Do not wring

MACHINE WASH

- Use a 40 wash cycle
- Use a non-biological mild detergent, **DO NOT** use fabric conditioner in any format
- Use a laundry bag if available



Hand Wash



Machine Wash at 40°C



Do Not Bleach



Iron Low Heat



Dry Clean



Tumble Dry Low Setting



Do Not Wring



Dry Flat



Drip Dry