



Sensory HUG

Jobskin®

Unit 13a Harrington Mill
Leopold Street, Long Eaton Nottingham, NG10 4QG UK
t: +44 (0) 115 973 4300 e: orders@jobskin.co.uk

www.jobskin.co.uk



EC REP

Advena Ltd. Tower Business Centre, 2nd Floor, Tower Street, Swatar, BKR 4013 Malta



Sensory Mini Hug

PCP33S

PCP33M

PCP33L

Sensory Maxi Hug

PCP34S

PCP34M

PCP34L

PCP34XL

The Sensory Hug can be used as part of a daily routine such as before or after school/work/daily activities, start of lessons, before going to sleep, or tied to specific activities such as eating, playing or going out. We suggest that you consult a suitably qualified person to recommend the optimal wearing regime based on the needs of the wearer.

Fitting Instructions:

1. Adjust any shoulder fastenings to fit the user.
2. Place the garment over the head with the fastenings at the front.
3. Fasten zip at required compression without zip strip, if necessary add 3cm, 5cm or 7cm or any combination up to an additional 15cm dependent on required compression.

User Instructions:

- Follow the above sequence for applying and removing the garment.
- Ensure that the fitting and effectiveness of this garment is periodically checked by a suitably qualified person.
- This garment can be worn under clothing or over clothing.
- Should any skin irritation occur when the garment is worn under clothing discontinue wearing and contact your treating clinician/provider.
- If the garment becomes tight or slack or damaged it will need to be replaced.
- As the name suggests the garment needs to provide the necessary compression to meet the needs of the user and may require some time to build up tolerance to wearing.

Care Instructions:

Wash your garment regularly, laundering will rejuvenate the material and help your garment to last longer. Use a non-biological mild detergent, do not use fabric softener. Do not dry in direct sunlight, on a radiator or in a tumble dryer. Dry naturally in a warm place. Wash dark colours separately.



Hand wash



Do not bleach



Dry flat



Do not tumble dry



Machine wash at 40°C



Do not iron



Drip dry



Do not wring