



Head Order Form

Order No.: _____ Patient Reference No.: _____



Circumference Measurements

(cm)

2	Around head at chin angle to crown of head	
5	At eyebrow level around head	
6	Around neck	



Length Measurements

Left (cm) Right (cm)

1	Corner of eye to corner of mouth		
3	Front of ear to ear around back of the head		
4	Distance from corner of eye to front of ear		
7	Bottom lip to point of chin (contour)		
8	Bottom chin to throat to neck finished length (contour)		
9	Broadest part bottom of nose (contour)		
10	Normal hairline to base of neck or finished length		

Please tick all that apply:

PG32 - Full Face Mask: Required Measurements: 1-8, 10

Closure: Zipper Velcro

Nose covering: Required Measurement: 9

Eyes to be cut at fitting

Ears to be marked for cutting at fitting

PG33 - Open Face Mask: Required Measurements: 2, 3, 5-8, 10

Closure: Zipper Velcro

Ears to be marked for cutting at fitting

PG34 - Wide Chin Strap: Required Measurements: 2, 5, 6, 7, 8

Velcro adjustable top of head (optional)

Ears to be marked for cutting at fitting

Specific Instructions