



Abdominal Support Order Form

All fields are required in order to process your order

Order Details

Date: _____
 Order No.: _____
 Contact Name: _____
 Contact Phone No.: _____
 Email: _____
 Hospital/Clinic: _____
 Delivery Address: _____

 _____ Post Code: _____

Please select:

Please note: we cannot process your order without indicating your options below.

Fabric:

White Coutil Beige Coutil White Brocade

Opening:

Right Side Left Side Centre Front

Stoma Aperture:

No Yes (Please mark position when fitting)

Standard features

The items below are standard with all Abdominal Supports, please indicate which options you require by ticking the relevant boxes.

Soft Fit: Yes

Plush Lining: Yes

Extra Large Fulcrum: Yes

Under Belt: Yes

Full Depth Elastic: Yes

Velcro D. Ring Fastening: Yes

Bones: Yes

Hernia Pad: Yes Extra Pad (chargeable)

Fat Restraining Pad: Yes Extra Pad (chargeable)

PLEASE NOTE:

- We only manufacture corsets which fall within the scope of this order form.
- Orders will be produced according to the supplied measures and features selected on the order form only. Corsets will then be supplied in a semi finished state for Soft-Fit.
- Please see the Corset adjustment guidelines on the reverse of the order form.

When completed, please click: customerservice@jobskin.co.uk to email your electronic order form

Please download your electronic forms directly from our website - www.jobskin.co.uk/file-download

Patient Details

Patient Reference No.: _____
 First Name: _____
 Surname: _____
 Date of Birth: _____

Circumference Measurements



F	Top edge	
G	Waist - in line with umbilicus	
H	Abdominal - in line with ASIS	
I	Hip - in line with trochanter	
J	Bottom - pubis to bottom of skirt	


Depth Measurements



A	Anterior depth - above waist	
B	Posterior depth - above waist	
C	Anterior depth - waist to pubis	
D	Posterior depth - waist to apex gluteus maximus	
E	Posterior depth - skirt	

Anterior

Posterior



Medigarments Corset Adjustment Guidelines and Prices

For Information only: this is a guideline for what charge may be incurred.

One Soft-Fit is included in the price of the corset.

Requested changes that fall within the conditions defined below are FOC following the first Soft-Fit.

- Aperture fitting
- Finishing of seams
- Shaping of bottom edge over hips
- Reduction in length of the velcro fastening straps

Some adjustments cannot be accommodated within the originally produced garment and will incur additional charges at soft fit.

Medium adjustments chargeable post Soft-Fit

- Reduction of circumferences up to 5cm
- Reduction of height above or below waist, centre back panel only - up to 5cm in total
- Increase of aperture, making larger only. Not a change of location

Major adjustments chargeable post Soft-Fit

- Increase of circumference up to 5cm, to be added to the centre back panel only
- Reduction of height above or below the waist, to include centre back, sides and centre front - up to 5cm in total

Adjustments which will result in a full re-cut

- Increase of circumference over 5cm
- Reduction of circumference over 5cm
- Fastening location change
- Aperture location change
- Increase of height above or below waist
- Reduction of height above or below the waist, to include centre back, sides and centre front - over 5cm

If requested, we will do a second Soft-Fit following the initial adjustments, however a second Soft-Fit will always be chargeable.

Corsets supplied finished after Soft-Fit adjustments will no longer be alterable

Medigarments non-conformity

If the garment has not been manufactured as per the original order instructions (a non-conformity), Medigarments will alter / replace the garment free of charge.