



Abdominal Support Order Form

All fields are required in order to process your order

Order Details

Date: \_\_\_\_\_  
Order No.: \_\_\_\_\_  
Contact Name: \_\_\_\_\_  
Contact Phone No.: \_\_\_\_\_  
Email: \_\_\_\_\_  
Hospital/Clinic: \_\_\_\_\_  
Delivery Address: \_\_\_\_\_  
Post Code: \_\_\_\_\_

Please select:

Please note: we cannot process your order without indicating your options below.

Fabric:

☐ White Coutil ☐ Beige Coutil ☐ White Brocade

Opening:

☐ Right Side ☐ Left Side ☐ Centre Front

Stoma Aperture:

☐ No ☐ Yes (Please mark position when fitting)

Standard features

The items below are standard with all Abdominal Supports, please indicate which options you require by ticking the relevant boxes.

Soft Fit: ☐ Yes  
Plush Lining: ☐ Yes  
Extra Large Fulcrum: ☐ Yes  
Under Belt: ☐ Yes  
Full Depth Elastic: ☐ Yes  
Velcro D. Ring Fastening: ☐ Yes  
Bones: ☐ Yes  
Hernia Pad: ☐ Yes ☐ Extra Pad (chargeable)  
Fat Restraining Pad: ☐ Yes ☐ Extra Pad (chargeable)

PLEASE NOTE:

- We only manufacture corsets which fall within the scope of this order form.
  - Orders will be produced according to the supplied measures and features selected on the order form only. Corsets will then be supplied in a semi finished state for Soft-Fit.
  - Please see the Corset adjustment guidelines on the reverse of the order form.
- When completed, please click: [customerservice@jobskin.co.uk](mailto:customerservice@jobskin.co.uk) to email your electronic order form

Please download your electronic forms directly from our website - [www.jobskin.co.uk/file-download](http://www.jobskin.co.uk/file-download)

Patient Details

Patient Reference No.: \_\_\_\_\_  
First Name: \_\_\_\_\_  
Surname: \_\_\_\_\_  
Date of Birth: \_\_\_\_\_

Circumference Measurements



F	Top edge	
G	Waist - in line with umbilicus	
H	Abdominal - in line with ASIS	
I	Hip - in line with trochanter	
J	Bottom - pubis to bottom of skirt	


Please Note: The maximum circumference we can manufacture is 150cms (for any circumference F, G, H, I, or J)

Depth Measurements



A	Anterior depth - above waist	
B	Posterior depth - above waist	
C	Anterior depth - waist to pubis	
D	Posterior depth - waist to apex gluteus maximus	
E	Posterior depth - skirt	

Please Note: The maximum total length we can manufacture is 60cms (lengths Anterior = A+C or Posterior = B+D+E)



Please Note: The maximum circumference we can manufacture is 150cms (for any circumference F, G, H, I, or J). The maximum total length we can manufacture is 60cms (lengths Anterior = A+C or Posterior = B+D+E)

## Medigarments Corset Adjustment Guidelines and Prices

**For Information only: this is a guideline for what charge may be incurred.**

**One Soft-Fit is included in the price of the corset.**

**Requested changes that fall within the conditions defined below are FOC following the first Soft-Fit.**

- Aperture fitting
- Finishing of seams
- Shaping of bottom edge over hips
- Reduction in length of the velcro fastening straps

**Some adjustments cannot be accommodated within the originally produced garment and will incur additional charges at soft fit.**

### **Medium adjustments chargeable post Soft-Fit**

- Reduction of circumferences up to 5cm
- Reduction of height above or below waist, centre back panel only - up to 5cm in total
- Increase of aperture, making larger only. Not a change of location

### **Major adjustments chargeable post Soft-Fit**

- Increase of circumference up to 5cm, to be added to the centre back panel only
- Reduction of height above or below the waist, to include centre back, sides and centre front - up to 5cm in total

### **Adjustments which will result in a full re-cut**

- Increase of circumference over 5cm
- Reduction of circumference over 5cm
- Fastening location change
- Aperture location change
- Increase of height above or below waist
- Reduction of height above or below the waist, to include centre back, sides and centre front - over 5cm

**If requested, we will do a second Soft-Fit following the initial adjustments, however a second Soft-Fit will always be chargeable.**

**Corsets supplied finished after Soft-Fit adjustments will no longer be alterable**

### **Medigarments non-conformity**

If the garment has not been manufactured as per the original order instructions (a non- conformity), Medigarments will alter / replace the garment free of charge.