

## Company Profile

Medigarments Ltd is the company behind the internationally recognised brands Jobskin® and SDO®. We lead the way in rehabilitation therapy and are the largest UK manufacturer of specialist custom compression garments supplied through therapy and orthotic professionals. Our reputation for the design, and supply of technically advanced garments has been built over 40 years supplying the highest quality made-to-measure products.

Our garments are used to manage scarring as a result of burns or trauma injury. In the field of dynamic compression, they are used to increase sensory and proprioceptive feedback as well as provide musculoskeletal support.

We are dedicated to providing the ultimate in therapy management, with the focus on improving patient quality of life.

## Medigarments Ltd

DESIGNED AROUND YOU

Medigarments Ltd  
Unit 13a Harrington Mill,  
Leopold Street,  
Long Eaton,  
Nottingham NG10 4QG  
United Kingdom

t: +44 (0) 115 973 4300  
e: customerservice@jobskin.co.uk

[www.jobskin.co.uk](http://www.jobskin.co.uk)

Please contact customer service for more details and information regarding our training courses.

V1 09/2022 | Jobskin® is a registered trade mark of Medigarments Ltd



Designed & Manufactured  
in the UK



# Designed Around You

Giving the best possible  
position for rehabilitation

## Jobskin® Gaiters



# Jobskin® Gaiters

Our range of made-to-measure and ready-to-wear gaiters deliver a high standard of craftsmanship in a class 1 medical device. Designed and manufactured using traditional techniques in our UK based factory, we are renowned for our quality and workmanship, where each product is crafted with passion, care, and the highest attention to detail.



## Product Description

Gaiters are designed as an external orthosis with structural inserts, which are used to immobilise the affected joint and hold the limb in a controlled extension or best possible position for rehabilitation. They provide a prolonged stretch to help maintain the gains achieved by manipulation, corrective surgery, reconstructive procedures or contractures. Designed to meet the individual's functional needs rather than to treat an individual condition or pathology, gaiters can be used by infants through to adulthood and are prescribed for use on one or more limbs.

## Indications

- Post injury / surgery
- Post fracture and removal of the cast
- Strains, sprains, and dislocations
- Juvenile Rheumatoid Arthritis
- Contracture management
- Immobilisation
- Spasticity management
- Post-operative corrective surgery or reconstruction



## Features

- Flexion gaiters made to customer requirements
- Angles of 5° increments from 5° to 90°
- Individual specified lengths
- Deluxe or standard linings for comfort and fit
- Available in several fabric options
- Simple application with fully adjustable straps