



All fields are required in order to process your order

Order Details

Date: _____ Order No.: _____
 Contact Name: _____
 Contact Phone No.: _____
 Email: _____
 Hospital/Clinic: _____
 Delivery Address: _____

 Post Code: _____

Patient Details

Patient Reference No.: _____
 First Name: _____
 Surname: _____
 Year of Birth: _____
 Please indicate: ☐ Male ☐ Female
 Please indicate: ☐ New Patient ☐ Existing Patient
 Diagnosis: _____

When all applicable fields are completed, please click:
customerservice@jobskin.co.uk to email your electronic order form

Please download your electronic forms directly from our website - www.jobskin.co.uk/file-download

All fabric choices are now offered as predefined themes, with colour-coordinated zips and threads. Please select one of the themes listed below:

Premium Original

Plain Powernet:

☐ Beige ☐ Tan ☐ Blossom ☐ Red ☐ Raspberry
☐ Classy Blue ☐ Denim Blue ☐ Black

Printed Powernet:

☐ Unicorn ☐ Safari Car ☐ Paw Print ☐ Pink Camo ☐ Green Camo
☐ Blue Camo ☐ Rainbow Unicorn ☐ Racing Car ☐ Mermaid ☐ 3D Football
☐ Green Football ☐ Flower Garden ☐ Picnic ☐ Abstract ☐ Pixel Brown
☐ Pixel Green

Premium Active - 50 UPF

☐ Eucalyptus Green ☐ Black

Premium Q10 - Q10 cosmetic ingredient

Plain Q10:

☐ Type 2 (white, fair) ☐ Type 3 (medium white to olive) ☐ Type 4 (olive, moderate brown) ☐ Type 5 (brown to dark brown) ☐ Type 6 (brown, very dark, brown to black)

Printed Q10:


☐ Fairy & Castle ☐ Dinosaurs

Order No.: _____ Patient Reference No.: _____


First select the required garment, then tick the box to confirm the standard configuration, or choose from the style options on page 5:

Garment (please indicate)

☐ PO 0535
Glove up to 7.6cm above wrist
☐ Left ☐ Right



☐ PO 0534
Glove above 7.6cm, to axilla or required finished length
☐ Left ☐ Right




Standard configurations for Gloves

Feature	☐ Gloves up to 14cm palm circumference (0535)	☐ Gloves above 14cm palm circumference (0535, 0534)
Proximal	Elastic cuff 1.5cm	Elastic cuff 2.5cm
Fingertips	Open / overlock	Open / overlock
Zipper	None / pull-on	None / pull-on


The following garments require style options from page 5:

Garment (please indicate)

☐ PO 0025
Individual Finger Sleeve Glove
(One finger only)
☐ Left ☐ Right



☐ PO 0036
Interdigital Web Spacer Glove - worn over a glove
☐ Left ☐ Right




Please note: you may use one form to order both left and right gloves, however they must be the same garment type and theme/fabric

Glove Order Form


Order No.: _____ Patient Reference No.: _____

Please use this outline in conjunction with the guide and grids on page 5

Use a traditional tape measure to record measurements. Type them in the boxes below or in the corresponding grid on the next page, both fill simultaneously.




Left Glove



Dorsal = Back of hand/arm
Volar = Palm side of hand/arm

Left Arm For Glove above wrist




Glove Order Form


Order No.: _____ Patient Reference No.: _____

Please use this outline in conjunction with the guide and grids on page 5


Use a traditional tape measure to record measurements. Type them in the boxes below or in the corresponding grid on the next page, both fill simultaneously.



Right Glove



Right Arm For Glove above wrist



Glove Order Form

Order No.: _____ Patient Reference No.: _____

Please note: Glove with less than 5 digits:
If you do not require all 5 digits to be included, please indicate by placing an 'X' in the measurement box.

Circumference Measurements		Left (cm)	Right (cm)
1	Thumb IP joint		
3	Index finger DIP		
4	Index finger PIP		
6	Middle finger DIP		
7	Middle finger PIP		
9	Ring finger DIP		
10	Ring finger PIP		
12	Little finger DIP		
13	Little finger PIP		
9a	Wrist crease		
22	Palm (ensure palm is opened out fully)		
26	Forearm at finished length of glove		

Length Measurements		Left (cm)	Right (cm)
10	From elbow joint crease to axilla or required length		
11	From wrist to elbow joint crease or required length		
15	Thumb web to finished length		
16	Index finger		
17	Middle finger		
18	Ring finger		
19	Little finger		
20	Thumb web to base of index finger		
21	Wrist to 2nd web		
21a	Wrist to 3rd web		
23	Wrist to 4th web		
25	Wrist to finished length		

For PO 0534 Glove above 7.6cm to axilla or required finished length

Limb Circumference Measurements		Left (cm)	Right (cm)
4	Top of arm level with axilla		
5	Mid upper arm		
6	Elbow joint crease, arm extended		
7	Upper forearm		
8	Lower forearm		
9a	Wrist crease		

PO 0025 - Individual Finger Garment

Required measurements:

Circumferences: 9a, 22

Lengths: 20, 21, 21a, 23

Please note: Circumferences and finished length of the required digit only

PO 0036 - Interdigital Web Spacer Garment

Required measurements:

Circumferences: 1, 4, 7, 9a, 10, 13, 22

Lengths: 15, 20, 21, 21a, 23

Style options when standard configuration is not selected

		Left	Right
Proximal elastic:	Overlock (no elastic)	☐	☐
	Regular (inverted) 1.5cm	☐	☐
	Regular (inverted) 2.5cm	☐	☐
	Cuff 1.5cm	☐	☐
	Cuff 2.5cm	☐	☐
	Silicone Regular (inverted) 2.5cm	☐	☐
	Silicone cuff 2.5cm	☐	☐
	Regular (inverted) 5cm (above elbow only)	☐	☐
	Cuff 5cm (above elbow only)	☐	☐
	Silicone Regular (inverted) 5cm (above elbow only)	☐	☐
	Silicone cuff 5cm (above elbow only)	☐	☐

Closed tips (tick if required, standard is open)		Left	Right
15	Thumb	☐	☐
16	Index finger	☐	☐
17	Middle finger	☐	☐
18	Ring finger	☐	☐
19	Little finger	☐	☐

Modifications available for all garments

All the following items will be an additional charge

Zippers - 1145 (tick if required)

Glove to wrist			
Zipper placement	☐ Inside of fabric	☐ Outside of fabric	
Position (please select: Ulnar, Radial, Dorsal, Volar)	Left	Right	
	☐	☐	
See measurement diagrams for guidance			
Length:	_____ cm		
Hook & eye (on fly behind the zip to assist donning)	☐	☐	

Inset Zippers - 1144 (tick if required)

Use placement pad to mark position if required			
Arm			
Extremity	☐ Distal	☐ Proximal	
Zipper placement	☐ Inside of fabric	☐ Outside of fabric	
Position (please select: Ulnar, Radial, Dorsal, Volar)	Left	Right	
	☐	☐	
See measurement diagrams for guidance			
Length:	_____ cm		

Order No.: _____ Patient Reference No.: _____

Sleeve Linings

Description	Product Code	Left	Right
Inner elbow lining (to protect fragile skin and provide comfort if required)	1167	☐	☐
Full elbow lining (as above)	1168	☐	☐

Reinforcements

Description	Product Code	Left	Right
Reinforced palm (outside of glove for high wear area)	0021	☐	☐
Non-slip grip to palm (glove only)	1197	☐	☐

Slant Inserts

Description	Product Code	Left	Right
Slant Inserts (a seam is sewn between the digits when additional pressure is required into the web spaces of the hand)	1169	☐	☐

Silon-TEX® II Insert

Silon-TEX® II fabric (sewn into garment)	1191	☐	Use placement pad to mark position
--	------	--------------------------	------------------------------------

Pockets & Pads

Pocket (sewn in for the insertion of pads to apply extra pressure to certain areas) Please specify position	0027	☐	Use placement pad to mark position
Silon-TEX® II pocket (as above but covered with Silon-TEX® II fabric)	1147	☐	Use placement pad to mark position

Foam Pads (to insert into pocket, please select foam thickness)

Low profile 5mm	1178	☐
Low density 20mm	1179	☐