



Gauntlet Order Form

All fields are required in order to process your order

Order Details

Date: _____ Order No.: _____
 Contact Name: _____
 Contact Phone No.: _____
 Email: _____
 Hospital/Clinic: _____
 Delivery Address: _____

 Post Code: _____

Patient Details

Patient Reference No.: _____
 First Name: _____
 Surname: _____
 Year of Birth: _____
 Please indicate: ☐ Male ☐ Female
 Please indicate: ☐ New Patient ☐ Existing Patient
 Diagnosis: _____

When all applicable fields are completed, please click:
customerservice@jobskin.co.uk to email your electronic order form

Please download your electronic forms directly from our website - www.jobskin.co.uk/file-download

All fabric choices are now offered as predefined themes, with colour-coordinated zips and threads. Please select one of the themes listed below:

Premium Original

Plain Powernet:

☐ Beige ☐ Tan ☐ Blossom ☐ Red ☐ Raspberry
☐ Classy Blue ☐ Denim Blue ☐ Black

Printed Powernet:

☐ Unicorn ☐ Safari Car ☐ Paw Print ☐ Pink Camo ☐ Green Camo
☐ Blue Camo ☐ Rainbow Unicorn ☐ Racing Car ☐ Mermaid ☐ 3D Football
☐ Green Football ☐ Flower Garden ☐ Picnic ☐ Abstract ☐ Pixel Brown
☐ Pixel Green

Premium Active - 50 UPF

☐ Eucalyptus Green ☐ Black

Premium Q10 - Q10 cosmetic ingredient

Plain Q10:

☐ Type 2 (white, fair) ☐ Type 3 (medium white to olive) ☐ Type 4 (olive, moderate brown) ☐ Type 5 (brown to dark brown) ☐ Type 6 (brown, very dark, brown to black)

Printed Q10:



☐ Fairy & Castle ☐ Dinosaurs

Gauntlet Order Form

Order No.: _____ Patient Reference No.: _____

First select the required garment, then tick the box to confirm the standard configuration, or choose from the style options on page 5:

Garment (please indicate)

☐ PO 0505 Gauntlet up to 7.6cm above wrist ☐ Left ☐ Right 	☐ PO 0536 Gauntlet above 7.6cm, to axilla or required finished length ☐ Left ☐ Right 
---	---

Standard configurations for Gauntlet

Feature	☐ Gauntlet up to 14cm palm circumference (0505)	☐ Gauntlet above 14cm palm circumference (0505, 0536)
Proximal	Elastic cuff 1.5cm	Elastic cuff 2.5cm
Distal	Overlock	Overlock
Zipper	None / pull-on	None / pull-on


Please note: you may use one form to order both left and right gauntlets, however they must be the same garment type and theme/fabric

Gauntlet Order Form

Order No.: _____ Patient Reference No.: _____


Please use this outline in conjunction with the guide and grids on page 5

Use a traditional tape measure to record measurements. Type them in the boxes below or in the corresponding grid on the next page, both fill simultaneously.




Left Gauntlet

Gauntlets are supplied with an open thumb tip



Dorsal = Back of hand/arm
Volar = Palm side of hand/arm

Left Arm For Gauntlet above wrist




Gauntlet Order Form

Order No.: _____ Patient Reference No.: _____


Please use this outline in conjunction with the guide and grids on page 5

Use a traditional tape measure to record measurements. Type them in the boxes below or in the corresponding grid on the next page, both fill simultaneously.




Right Gauntlet

Gauntlets are supplied with an open thumb tip



Right Arm
For Gauntlet
above wrist



Gauntlet Order Form

Order No.: _____ Patient Reference No.: _____

Circumference Measurements

		Left (cm)	Right (cm)
1	Thumb IP joint		
9	Wrist crease		
22	Palm (ensure palm is opened out fully)		
26	Finished length		

Length Measurements

		Left (cm)	Right (cm)
10	From elbow joint crease to axilla or required length		
11	From wrist to elbow joint crease or required length		
15	Thumb web to finished length		
20	* Thumb web to end of gauntlet		
21	* Wrist to end of gauntlet		
23	* Wrist to end of gauntlet		
25	Wrist to finished length of gauntlet		

* Please Note: The lengths required for this style need to be to the length required and NOT to the web spaces.

Limb Circumference Measurements

		Left (cm)	Right (cm)
4	Top of arm level with axilla		
5	Mid upper arm		
6	Elbow joint crease, arm extended		
7	Upper forearm		
8	Lower forearm		
9	Wrist crease		

Style options when standard configuration is not selected

configuration is not selected		Left	Right
Proximal elastic:	Overlock (no elastic)	☐	☐
	Regular (inverted) 1.5cm	☐	☐
	Regular (inverted) 2.5cm	☐	☐
	Cuff 1.5cm	☐	☐
	Cuff 2.5cm	☐	☐
	Silicone Regular (inverted) 2.5cm	☐	☐
	Silicone cuff 2.5cm	☐	☐
	Regular (inverted) 5cm (above elbow only)	☐	☐
	Cuff 5cm (above elbow only)	☐	☐
Silicone Regular (inverted) 5cm (above elbow only)	☐	☐	
	Silicone cuff 5cm (above elbow only)	☐	☐
Distal elastic:	Overlock (no elastic)	☐	☐
	Regular (inverted) 1.5cm	☐	☐
	Cuff 1.5cm	☐	☐

Modifications available for all garments

All the following items will be an additional charge

Zippers - 1145 (tick if required)

Gauntlet to wrist		
Zipper placement	☐ Inside of fabric	☐ Outside of fabric
Position (please select: Ulnar, Radial, Dorsal, Volar)	Left ☐	Right ☐
See measurement diagrams for guidance		
Length: _____ cm		
Hook & eye (on fly behind the zip to assist donning)	☐	☐

Inset Zippers - 1144 (tick if required)

Use placement pad to mark position if required			
Arm	☐	☐	
Extremity	☐ Distal	☐ Proximal	
Zipper placement	☐ Inside of fabric	☐ Outside of fabric	
Position (please select: Ulnar, Radial, Dorsal, Volar)	Left ☐	Right ☐	
See measurement diagrams for guidance			
Length: _____ cm			

Sleeve Linings

Description	Product Code	Left	Right
Inner elbow lining (to protect fragile skin and provide comfort if required)	1167	☐	☐
Full elbow lining (as above)	1168	☐	☐

Reinforcements

Description	Product Code	Left	Right
Non-slip grip to palm (glove only)	1197	☐	☐

Silon-TEX® II Insert

Silon-TEX® II fabric (sewn into garment)	1191	☐ Use placement pad to mark position
--	------	---

Pockets & Pads

Pocket (sewn in for the insertion of pads to apply extra pressure to certain areas) Please specify position	0027	☐ Use placement pad to mark position
Silon-TEX® II pocket (as above but covered with Silon-TEX® II fabric)	1147	☐ Use placement pad to mark position

Foam Pads (to insert into pocket, please select foam thickness)

Low profile 5mm	1178	☐
Low density 20mm	1179	☐