



Arm Gaiter Order Form  
Made-to-measure

All fields are required in order to process your order

Order Details

Date: \_\_\_\_\_  
Order No.: \_\_\_\_\_  
Contact Name: \_\_\_\_\_  
Contact Phone No.: \_\_\_\_\_  
Email: \_\_\_\_\_  
Hospital/Clinic: \_\_\_\_\_  
Delivery Address: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
Post Code: \_\_\_\_\_

Patient Details

Patient Reference No.: \_\_\_\_\_  
First Name: \_\_\_\_\_  
Surname: \_\_\_\_\_  
Date of Birth: \_\_\_\_\_

PLEASE NOTE: Gaiters are non-returnable.  
Please ensure accurate measurements as we  
are unable to make alterations.

One Garment Per Form

Specifications

Arm: ☐ Left ☐ Right  
Style: ☐ KNS08 Extension ☐ KNS08F Flexion\*  
\* Please indicate degree of flexion required: \_\_\_\_\_  
Flexion is available in increments of 5 degrees, beginning at zero.  
Incorrectly supplied flexion angles will automatically be round-  
ed up or down to the nearest 5 degrees to a maximum of 90  
degrees.  
Straight steels: ☐ 2 ☐ 3 ☐ 4  
Please note: Flexion gaiters are supplied with 2 flexion steels

Circumference  
Measurements



A	Proximal circumference	
B	Distal circumference	

Length  
Measurements




C	Proximal end to elbow joint (Measure Posterior surface)	
D	Distal end to elbow joint (Measure Posterior surface)	
G	Total length of extension gaiter (Measure Medial Border)	

Fabrics

☐ Bumblebee ☐ Camouflage ☐ Dinosaurs  
☐ Hearts & Flowers ☐ Black ☐ Rainbow Unicorns  
☐ Urban Camo ☐ Denim ☐ Claret

Lining: ☐ Cotton ☐ Towelling & Foam Padded (chargeable)



When completed, please click: [customerservice@jobskin.co.uk](mailto:customerservice@jobskin.co.uk) to email your electronic order form

Please download your electronic forms directly from our website - [www.jobskin.co.uk/file-download](http://www.jobskin.co.uk/file-download)