




Classic

MTM Pressure Garment



Pressure Garments

Made to Measure order forms

Order form: Foot

Guidelines for Foot order form

Recording measurements

- **Please record your measurements using the outline drawing, measurement boxes or both**
- All measurements are recorded in metric. For example 12.4cm

Sock

- Zippers can be inserted if the garment is difficult to apply
- Please consider the finished length of your garment in order to minimise the garment from slipping and to ensure comfort for your patient
- Sock can be closed with powernet fabric or soft cotton
- Open toe sock: it is advisable to finish at the metatarsal heads
- For scarring to the toes or into the webs, a foot glove should be considered
- Please contact our therapy specialist or the manufacturing unit for assistance with measuring

Silicone textile inserts or pad inserts

- Please use the body outline tracing pad provided to indicate any special requirements or to mark the area where you would like a silicone textile insert or a pad inserted
- Please mark clearly using a black pen
- Alternatively these can also be marked at the 1st fitting for accuracy if desired


Patient Name:
Date:
Product Name:
Quantity:



Jobskin

Classic

MTM Pressure Garment



	Measurements	Left (cm)	Right (cm)
15:	Knee joint		
16:	Upper calf		
17:	Mid calf		
18:	Lower calf		
19:	Upper margin of medial malleolus		
21:	Medial malleolus to required top of sock		
24:	Around foot and heel under malleolus		
25:	Instep or waist of foot		
26:	Foot at metatarsal heads		
27:	Upper margin medial malleolus to sole foot		
28:	Metatarsal heads to heel (medial)		
29:	Metatarsal heads to tip of toes (medial)		

Sock style/options (tick if required)	L	R
Open toe		
Enclosed toe (please circle): self - soft		
Ankle contracture seam		
Reinforced heel		
Zipper location (please circle): lateral - medial		

Specific instructions:
