


PATIENT REFERENCE NUMBER OR NAME:

Product Name:

Quantity:

Date:

Jobskin Ltd: Unit 13a Harrington Mill, Leopold Street, Long Eaton, Nottingham, NG10 4QG.
t: +44 (0) 115 973 4300 | f: +44 (0) 115 973 3902 | e: orders@jobskin.co.uk | w: www.jobskin.co.uk



Measurement for Circumferences		Left (cm)	Right (cm)
1	Thumb IP joint		
3	Index finger DIP		
4	Index finger PIP		
6	Middle finger DIP		
7	Middle finger PIP		
9	Ring finger DIP		
10	Ring finger PIP		
12	Little finger DIP		
13	Little finger PIP		
22	Palm (ensure that the palm is opened out fully)		
24	Wrist crease		
26	Forearm at finished length of glove		

Measurement for Lengths		Left (cm)	Right (cm)
15	Thumb		
16	Index finger		
17	Middle finger		
18	Ring finger		
19	Little finger		
20	Thumb web to base of index finger		
21	Wrist to 2nd web		
21a	Wrist to 3rd web		
23	Wrist to 4th web		
25	Wrist to finished length of glove		

Measurements Required for Gauntlet Style
Circumferences: 1, 15, 22, 24, 25, 26
Lengths: 20, 21, 23, 25
Please Note: The lengths required for this style need to be to the length required and NOT the web spaces.

Close Tip (tick if required)	Left	Right
Thumb		
Index Finger		
Middle Finger		
Ring Finger		
Little Finger		

Web Spacers (please indicate)
Th _____ 2nd _____ 3rd _____ 4th _____
Spacers in Elastic (please indicate)
<input type="checkbox"/> 4mm width only
Spacers in Cotton Webbing (please indicate width required)
<input type="checkbox"/> 4mm <input type="checkbox"/> 8mm

Specific Instructions:

Glove/Gauntlet Style Options (tick if required)	Left	Right
Proximal elastic (please tick): <input type="checkbox"/> Regular <input type="checkbox"/> Cuff		
Proximal elastic, non-slip		
Reinforced palm (outside)		
Non-slip grip to palm		
Lined palm inside		
Velcro tab closure (proximal end)		
Popper closure (proximal end)		
Pad (please give details and enclose drawing)		
Zipper Location (please tick): <input type="checkbox"/> Ulnar <input type="checkbox"/> Radial <input type="checkbox"/> Dorsal		