


**Jobskin Ltd:** Unit 13a Harrington Mill, Leopold Street, Long Eaton, Nottingham, NG10 4QG, United Kingdom  
t: +44 (0) 115 973 4300 | f: +44 (0) 115 973 3902 | e: orders@jobskin.co.uk | w: www.jobskin.co.uk

<b>Product: PG32 - Full Face Mask</b> <i>(please tick)</i>	
<b>Measurements:</b> Nos. 1 - 8 and 10 <b>required</b> only	
<b>Closure</b> <i>(please tick)</i> : <input type="checkbox"/> Zipper <input type="checkbox"/> Velcro	
<b>Nose covering:</b> No. 9 required <i>(please tick)</i>	
Eyes to be cut at fitting	
Ears to be marked for cutting at fitting	

<b>Product: PG33 - Open Face Mask</b> <i>(please tick)</i>	
<b>Measurements:</b> Nos. 2, 3, 5-8 and 10 <b>required</b> only	
<b>Closure</b> <i>(please tick)</i> : <input type="checkbox"/> Zipper <input type="checkbox"/> Velcro	
Ears to be marked for cutting at fitting	

<b>Product: PG34 - Wide Chin Strap</b> <i>(please tick)</i>	
<b>Measurements:</b> Nos. 2, 5, 6, 7, 8 <b>required</b> only	
<b>Velcro adjustable top of head</b> <i>(please tick)</i>	
Ears to be marked for cutting at fitting	

<b>Product: PG31 - Chin Strap</b> <i>(please tick)</i>	
<b>Measurements:</b> Nos. 2, 5, 6, 7, 8 <b>required</b> only	
<b>Velcro adjustable top of head</b> <i>(please tick)</i>	



**PATIENT REFERENCE NUMBER OR NAME:**

**Product:** *(Please indicate in the product box which product you require)*

**Measurements:** *(Please indicate in the measurement boxes below)*

**Quantity:** **Date:**

Measurement for Circumferences		cm
2	Around head at chin angle to crown of head	
5	At eyebrow level around head	
6	Around neck	
Measurement for Lengths		cm
1	Corner of eye to corner of mouth	
3	Front of ear to ear around back of the head	
4	Distance from corner of eye to front of ear	
7	Bottom lip to point of chin (contour)	
8	Bottom chin to throat to neck finished length (contour)	
9	Broadest part bottom of nose (contour)	
10	Normal hairline to base of neck/finished length	

<b>Specific Instructions:</b>	
-------------------------------	--