

# Guidelines for Lower Extremity order form

## Recording measurements

- Please record your measurements using the outline drawing, measurement boxes or both
- All measurements are recorded in metric. For example 12.4cm

## Pants/Shorts/Leg sleeve/One leg hip style garment order form

- Please select your design and then mark **only the required** measurements for each style of garment
- Measure your patient standing if possible
- Choose your style carefully according to your assessment and your patients lifestyle

## Waist

- The waist is where the patient wears their trousers, briefs or skirt and you will need to take your waist measurement at this point

## Fuller figure

- In order to obtain a good fit for the obese patient or a child with a protruding abdomen the pants/shorts may need to be higher at the front than the back
- Please consider this when taking your measurements and indicate this on your order form

## Children in nappies

- Please measure over the nappy

## Leg sleeve

- You will need to consider the finished length carefully in order to minimise the garment from slipping down and to ensure comfort for your patient
- You may also wish to consider the one leg hip style garment (PG 27) or the one leg chap style (PG 26) as an alternative style

## Silicone textile insert or pad insert

- Please use the body outline tracing pad provided or reverse of pad to indicate any special requirements or to mark the area where you would like a silicone textile insert or a pad inserted
- Please mark clearly using a black pen
- Alternatively these can also be marked at the first fitting for accuracy if desired

Premium MTM Data Order Forms can also be downloaded from our website [www.jobskin.co.uk](http://www.jobskin.co.uk)


<b>Patient Name:</b>
<b>Date:</b>
<b>Product Name:</b>
<b>Quantity:</b>

<b>PG57 and PG 58: Panty, chest height with straps additional measurements required</b>	<b>(cm)</b>
<b>B:</b> Chest at axilla level	
<b>C:</b> Waist to front height of garment <b>excluding scoop under arm</b>	
Length of straps (please state length/cm)	
Scoop under arm (standard depth 3cm)	<b>Yes / No</b>

<b>PG59: French Knicker Panty</b>	<b>(cm)</b>
<b>A:</b> Waist	
<b>J:</b> Hips	
<b>D:</b> Waist to hips	
Length above waist if required (High waisted only)	
Circumference above waist if required (high waisted only)	
<b>K:</b> Waist to waist through crotch	

<b>New Lace Option:</b> Please tick box if you require the new lace trim on the garment <input type="checkbox"/>
---

<b>Specific instructions:</b>



# Jobskin

	Measurements	Left (cm)	Right (cm)
<b>A:</b>	Waist		
<b>D:</b>	Waist to hips PG27 only		
<b>J:</b>	Hips		
<b>13:</b>	Top of thigh level with gluteal fold		
<b>13a:</b>	Upper thigh		
<b>13b:</b>	Mid thigh		
<b>14:</b>	Lower thigh		
<b>15:</b>	Knee joint (in line with mid patella)		
<b>16:</b>	Upper calf		
<b>17:</b>	Mid calf		
<b>18:</b>	Lower calf		
<b>19:</b>	Ankle at upper margin of medial malleolus		
<b>20:</b>	Inside leg (into groin) to knee joint or required length if above knee		
<b>21:</b>	Knee joint to ankle or required length		
<b>22:</b>	Waist to knee or required length of garment if above knee		

Style/options (tick if required)	L	R
Cotton lined crotch <b>or</b> open crotch (please circle)		
Velcro adjustable waist		
Double width waistband		
Distal end garment elastic cuff <b>or</b> overlock (please circle)		
Lining behind knee		
Reinforced knee		
New Boxer diagonal fly <b>or</b> zipper fly opening (please circle)		
Velcro tabs (to attach to vest) (please circle): YES		
<b>Extra options leg styles</b>		
Proximal elastic double width		
Proximal elastic non-slip		

Unit 13a Harrington Mill | Leopold Street | Long Eaton | Nottingham | NG10 4QG  
 T: +44 (0)115 973 4300 | F: +44 (0)115 973 3902 | W: www@jobskin.co.uk



PG 31 chin strap  
 PG 34 wide chin strap



PG 32 full mask



PG 33 open face mask



PG 13 glove to wrist

L  R



PG 16 gauntlet to wrist

L  R



PG 14 glove to elbow  
 L  R




PG 17 gauntlet above wrist  
 L  R



PG 10 single arm  
 L  R



PG 55 vest with straps



PG 56 vest one arm long or short sleeve



PG 57 panty chest height short leg  
 PG 58 panty long leg



PG 52 sternal vest



PG 1 vest no sleeves  
 PG 2 vest short long sleeves  
 Bra cups in Vest



PG 3 leotard no sleeves



PG 4 leotard short/long sleeves



PG 5 Single Elastic trunk belt  
 PG Double Elastic trunk belt



PG 12A arm sleeve shoulder flap  
 L  R



PG 59 French Knicker



PG 60 Pregnancy Panty



PG 20 open/closed toe sock  
 PG 25 leg sleeve any length  
 L  R



PG 21 below knee stocking  
 PG 22 thigh length stocking  
 L  R



PG 23 panty short leg  
 PG 24 panty long leg



PG 26 single leg with waist attachment  
 L  R



PG 27 single leg hip style  
 L  R

# Pressure Garments

## Made to Measure Order Forms

### Patient details

Surname:
First name:
Please indicate: Male <input type="checkbox"/> Female <input type="checkbox"/>
Date of birth:
Diagnosis:

### Customer details

Date:
Hospital and address:
Your Order No:
Therapist name:
Direct contact No:

**Powernet fabric:** White  Beige  Bronze  Pastel Pink  Pastel Blue  Purple  Black

Red  Royal Blue  Camouflage

**Polycotton fabric:** White  Beige  Pastel Pink  Purple  Navy Blue  Black

**Thread colour:** White  Beige  Bronze  Pastel Blue  Pastel Pink  Mint Green

Royal Blue  Denim Blue  Green  Red  Purple  Black  Bright Pink

**Motif number:** \_\_\_\_\_

**Delivery address:** \_\_\_\_\_

**Post code:** \_\_\_\_\_

**For express delivery, please contact customer service. Date appointed for \_\_/\_\_/20\_\_**

**When completed, please fax back the relevant pages (including this cover page) to 0115 973 3902**