

Guidelines for Lower Extremity order form

Recording measurements

- Please record your measurements using the outline drawing, measurement boxes or both
- All measurements are recorded in metric. For example 12.4cm

Pants/Shorts/Leg sleeve/One leg hip style garment order form

- Please select your design and then mark **only the required** measurements for each style of garment
- Measure your patient standing if possible
- Choose your style carefully according to your assessment and your patients lifestyle

Waist

- The waist is where the patient wears their trousers, briefs or skirt and you will need to take your waist measurement at this point

Fuller figure

- In order to obtain a good fit for the obese patient or a child with a protruding abdomen the pants/shorts may need to be higher at the front than the back
- Please consider this when taking your measurements and indicate this on your order form

Children in nappies

- Please measure over the nappy

Leg sleeve

- You will need to consider the finished length carefully in order to minimise the garment from slipping down and to ensure comfort for your patient
- You may also wish to consider the one leg hip style garment (PG 27) or the one leg chap style (PG 26) as an alternative style

Silicone textile insert or pad insert

- Please use the body outline tracing pad provided or reverse of pad to indicate any special requirements or to mark the area where you would like a silicone textile insert or a pad inserted
- Please mark clearly using a black pen
- Alternatively these can also be marked at the first fitting for accuracy if desired

Premium MTM Data Order Forms can also be downloaded from our website www.jobskin.co.uk


Patient Name:
Date:
Product Name:
Quantity:

PG57 and PG 58: Panty, chest height with straps additional measurements required	(cm)
B: Chest at axilla level	
C: Waist to front height of garment excluding scoop under arm	
Length of straps (please state length/cm)	
Scoop under arm (standard depth 3cm)	Yes / No

PG59: French Knicker Panty	(cm)
A: Waist	
J: Hips	
D: Waist to hips	
Length above waist if required (High waisted only)	
Circumference above waist if required (high waisted only)	
K: Waist to waist through crotch	

New Lace Option: Please tick box if you require the new lace trim on the garment <input type="checkbox"/>

Specific instructions:



Jobskin

	Measurements	Left (cm)	Right (cm)
A:	Waist		
D:	Waist to hips PG27 only		
J:	Hips		
13:	Top of thigh level with gluteal fold		
13a:	Upper thigh		
13b:	Mid thigh		
14:	Lower thigh		
15:	Knee joint (in line with mid patella)		
16:	Upper calf		
17:	Mid calf		
18:	Lower calf		
19:	Ankle at upper margin of medial malleolus		
20:	Inside leg (into groin) to knee joint or required length if above knee		
21:	Knee joint to ankle or required length		
22:	Waist to knee or required length of garment if above knee		

Style/options (tick if required)	L	R
Cotton lined crotch or open crotch (please circle)		
Velcro adjustable waist		
Double width waistband		
Distal end garment elastic cuff or overlock (please circle)		
Lining behind knee		
Reinforced knee		
New Boxer diagonal fly or zipper fly opening (please circle)		
Velcro tabs (to attach to vest) (please circle): YES		
Extra options leg styles		
Proximal elastic double width		
Proximal elastic non-slip		

Unit 13a Harrington Mill | Leopold Street | Long Eaton | Nottingham | NG10 4QG
 T: +44 (0)115 973 4300 | F: +44 (0)115 973 3902 | W: www@jobskin.co.uk



PG 31 chin strap
 PG 34 wide chin strap



PG 32 full mask



PG 33 open face mask



PG 13 glove to wrist

L R



PG 16 gauntlet to wrist

L R



PG 14 glove to elbow
 L R



PG 17 gauntlet above wrist
 L R



PG 10 single arm
 L R



PG 55 vest with straps



PG 56 vest one arm long or short sleeve



PG 57 panty chest height short leg
 PG 58 panty long leg



PG 52 sternal vest



PG 1 vest no sleeves
 PG 2 vest short long sleeves
 Bra cups in Vest



PG 3 leotard no sleeves



PG 4 leotard short/long sleeves



PG 5 Single Elastic trunk belt
 PG Double Elastic trunk belt



PG 12A arm sleeve shoulder flap
 L R



PG 59 French Knicker



PG 60 Pregnancy Panty



PG 20 open/closed toe sock
 PG 25 leg sleeve any length
 L R



PG 21 below knee stocking
 PG 22 thigh length stocking
 L R



PG 23 panty short leg
 PG 24 panty long leg



PG 26 single leg with waist attachment
 L R



PG 27 single leg hip style
 L R

Pressure Garments

Made to Measure Order Forms

Patient details

Surname:
First name:
Please indicate: Male <input type="checkbox"/> Female <input type="checkbox"/>
Date of birth:
Diagnosis:

Customer details

Date:
Hospital and address:
Your Order No:
Therapist name:
Direct contact No:

Powernet fabric: White Beige Bronze Pastel Pink Pastel Blue Purple Black

Red Royal Blue Camouflage

Polycotton fabric: White Beige Pastel Pink Purple Navy Blue Black

Thread colour: White Beige Bronze Pastel Blue Pastel Pink Mint Green

Royal Blue Denim Blue Green Red Purple Black Bright Pink

Motif number: _____

Delivery address: _____

Post code: _____

For express delivery, please contact customer service. Date appointed for __/__/20__

When completed, please fax back the relevant pages (including this cover page) to 0115 973 3902