

SDO[®] Garments for Inoperable Hernias

This leaflet is to assist Orthotists in the management of Inoperable Hernias where the standard binder or corset is not effective.

The style of garment required will be determined by the size and position of the hernia.

Please select one of the following garment styles:

- PCP29 - High waisted shorts with or without shoulder straps
- PCP01 - Sleeveless vest
- PCP03 - Sleeveless/short leg body suit

Additional information required with the garment order:

1. Photos of the hernia. These Images will be used to help with garment design and manufacture
2. Patient Consent Form, signed by the patient, for the above photos
3. Additional hernia measures: length, width, and depth of hernia lying in supine with as much reduction as possible
4. When completing circumferential body measures over the hernia take the measurements as firmly as possible to gain the best reduction
5. Request the bottom of vest to be sewn on stretch at the front of the garment (only for vest style option)
6. Dependent on size of hernia you may want to request an open side seam zip so that donning and doffing is made easier.

If you would like assistance or advice, please contact our
SDO Clinical Manager:

Lydia Dean

m: +44 (0) 7780 688020 | **e:** lydia.dean@jobskin.co.uk

Jobskin Ltd UK:

Unit 13a Harrington Mill, Leopold Street, Long Eaton, Nottingham, NG10 4QG, UK.
t: +44 (0) 115 973 4300 | **f:** +44 (0) 115 973 3902 | **e:** orders@jobskin.co.uk

Jobskin Ltd Ireland:

Unit C3 Innovation Works, Tipperary Technology Park, Thurles, Co Tipperary, Ireland.
t: +353 (0) 504 90100 | **f:** +353 (0) 504 22735 | **e:** ireorders@jobskin.ie

www.jobskin.co.uk



Designed & Manufactured
in the UK