

Shorts and Pants Order Form









Unit 13a Harrington Mill, Leopold Street, Long Eaton, Nottingham, NG10 4QG

t: +44 (0) 115 973 4300 | f: +44 (0) 115 973 3902 | e: orders@jobskin.co.uk | w: www.jobskin.co.uk

Please indicate specific garment by selecting the relevant tick box.

Type of SDO Required	Product No.	Quantity
Pants Long Leg: <input type="checkbox"/> Open Crotch <input type="checkbox"/> Closed Crotch	PCP15	
Pants Short Leg: <input type="checkbox"/> Open Crotch <input type="checkbox"/> Closed Crotch	PCP15A	
High Waisted Pants with Shoulder Straps	PCP29	
Measuring session fee	PCP20	

Reinforcement Panels for Shorts and Pants (Please select panel and tick box)

Standard Panels	Lumbar Panels	Hip Panels	Medial Upper Leg Panels
SP (Std Panels included with Garment) To assist hip external rotation	<input type="checkbox"/> LP To assist postural control	<input type="checkbox"/> HP To assist hip control and external rotation	<input type="checkbox"/> ULM To assist adduction and internal rotation
 <p>Anterior Posterior</p>	 <p>Anterior Posterior</p>	 <p>Anterior Posterior</p>	 <p>Anterior Posterior</p>
Anterior Pelvic Panels	Posterior Pelvic Panels	External Rotation Panels	Internal Rotation Panels
<input type="checkbox"/> APP To resist posterior tilt	<input type="checkbox"/> PPP To resist anterior tilt and encourage hip extension	<input type="checkbox"/> ERP To assist external rotation at the hips	<input type="checkbox"/> IRP To assist internal rotation at the hips
 <p>Anterior</p>	 <p>Posterior</p>	 <p>Anterior Posterior</p>	 <p>Anterior Posterior</p>

Please note: If you require any additional or bespoke panels to the above, please include in "Any other information" below, or attach separately. Please also see Panel Poster for a complete list of available panels on "Download" tab on www.jobskin.co.uk

Any other information, if required attach separately.

Delivery Address

Post code:

Patient Details
(Please complete Patient Ref. No. OR Name)

PATIENT REFERENCE NO.:

FIRST NAME:

SURNAME:

Date of Birth:

Please indicate: Male Female

New Patient Existing Patient

Diagnosis (must be completed):

Date:

Hospital/Clinic:

Purchasers Details

Address:

Post Code:

Consultant:

Name:

Measured by: OT PT Orthotist Jobskin

Hospital Order No.:

Contact Phone No.:

Personalisation Choices (Please indicate using code or name)

Fabric Colour:

Reinforcement Panel Colour:

Thread Colour:

Motif:

Binding:

Printed Zipper:

Wears AFO's: Yes No

Silicone Edging: Yes No

Zips for Pants:

Crotch: Short Leg only

Inside Leg: Short Leg only

Outside Leg: Short Leg only

Side Seam: Waist to Hip
 Hip to end of Leg

First SDO: Yes No

Fabric Colour		
White	Beige	Black
Pastel Pink	Navy Blue	

Reinforcement Panel Colour		
White	Beige	Black
Pastel Pink	Red	Bronze
Pastel Blue	Purple	Pink
Royal Blue		
Floral	Paisley	Camouflage
Marvel	Bull Dog	Unicorn

Thread Colour		
White	Beige	Black
Pastel Pink	Red	Bronze
Pastel Blue	Purple	Pink
Royal Blue	Yellow	Orange
Bright Pink	Navy	Denim
Mint Green	Green	Blue

When completed, please fax or email back to Jobskin - Fax: 0115 973 3902 | Email: orders@jobskin.co.uk

Shorts and Pants Order Form



(PLEASE COMPLETE PATIENT REFERENCE NO. OR NAME)

PATIENT REFERENCE NO.:

FIRST NAME: SURNAME:

Date of Assessment:

Please make sure when measuring for a SDO garment that the patient is measured lying down, supine and in the best position of alignment.



Measurements for Shorts and Pants - Circumference		Left (cm)	Right (cm)
A	Level with umbilicus 1cm less than "at rest" circumference.		
J	Level with Anterior Superior Iliac Spine (ASIS).		
13	Leg level with the groin and gluteal fold.		
13A	Mid-thigh level with muscle bulk 5mm less than "at rest" circumference.		
14	Above knee at suprapatellar margin or where shorts end.		
15	Level with knee joint (in line with mid patella).		
16	Below knee at upper margin of tibial tuberosity.		
17	Level with bulk of calf 5mm less than "at rest" circumference.		
18	Level with mid-lower leg, at point of musculotendinous junction.		
19	Level with upper margin of medial malleolus.		
Instep	Around instep or waist of the foot (only if stirrup required).		

Measurements for Shorts and Pants - Lengths		Left (cm)	Right (cm)
D	Point level with umbilicus to Anterior Superior Iliac Spine (ASIS).		
20	Inside leg groin to knee or end of SDO.		
21	Knee joint in line with mid-patella to upper margin of medial malleolus - measured on medial border.		
22	Point level with umbilicus on anterior surface to knee or end of shorts or long leg pants.		

Extra Measurements for High Waisted Shorts/Pants - Circumference		Left (cm)	Right (cm)
B	Top of shorts or pants.		

Extra Measurements for High Waisted Shorts/Pants - Lengths		Left (cm)	Right (cm)
C	Top of shorts to waist (A)		