


Design Choices

Nose Covering	<ul style="list-style-type: none"> - Provides pressure to bridge and tip of nose. - A contour splint is recommended under full face mask possibly with padding. - Give measurements A and B on the Jobskin Order Form 502.
Open Head Mask	<ul style="list-style-type: none"> - This style is open from the eyebrows to the same distance at the back of the head. - There is a strap from centre back to centre back.

Modifications



0014 Closed Ear Flap	<ul style="list-style-type: none"> - Can be incorporated into face mask or modified chin strap.
0017 Head Band	<ul style="list-style-type: none"> - Separate head band. - Measurement number 2 and width required.
1166 Lip Covering	<ul style="list-style-type: none"> - A detachable lip covering to provide pressure to the lips. - Can be removed for eating. - Can be included on both face mask and chin strap.
0034 Detachable Turtle Neck	<ul style="list-style-type: none"> - Made from elastic, 2.5cm (1"), 5cm (2") or 7.62cm (3") wide. - Used to provide pressure.
0027 Pockets	<ul style="list-style-type: none"> - This modification allows for the insertion of padding in areas where pressure is difficult to achieve. - The position and size must be indicated on face drawing form 502 or on a piece of paper